

BIOTOOLS

BIOTOOLS B&M LABS.S.A.

Manufactured by:

BIOTOOLS B&M Labs, S.A.
Valle de Tobalina - 52 - Nave 39
28021 Madrid
Spain

Tel. (34) 91 710 00 74
Fax (34) 91 505 31 18
e-mail: info@biotools.eu
www.biotools.eu

SPEEDTOOLS TISSUE DNA EXTRACTION KIT

*Kit for the Isolation of Genomic DNA from
Tissue, Cells, Bacteria, Yeast and Body Fluids*

Instructions for Use (Ref. 21.135M/6/7)

**PLEASE READ THE INSTRUCTIONS FOR USE THOROUGHLY BEFORE USING THE KIT,
ESPECIALLY IF YOU ARE NOT FAMILIAR WITH THE PROTOCOL.**

Table of Contents

Instructions for Use	1
2. KIT CONTENTS.....	2
3. INTENDED USE.....	2
4. PREPARATION OF WORKING SOLUTIONS AND STORAGE CONDITIONS	3
5. DNA ELUTION PROCEDURES	4
6. INSTRUCTION FOR USE.....	5
A. STANDARD PROTOCOL (for human or animal tissue and cultured cells).....	5
B. Protocol for mouse or rat tails	7
C. Protocol for bacteria	8
D. Protocol for yeast.....	9
E. Protocol for dried blood spots (e.g. Guthrie cards; NucleoCards; and FTA® cards).....	10
F. Protocol for genomic DNA and viral DNA from blood samples.....	11
G. Protocol for hair roots	12
H. Protocol from paraffin-embedded tissue	12
I. Protocol for detection of <i>Mycobacterium tuberculosis</i> or <i>Legionella pneumophila</i> in sputum or bronchoalveolar lavage.....	13
J. Protocol for extraction of genomic DNA from stool	13
K. Protocol for viral DNA from stool.....	14
L. Protocol for detection of EHEC bacteria in food (e.g. fresh cows' milk)	15
M. Protocol for purification of bacterial DNA (e.g. <i>Borrelia burgdorferi</i>) from urine.....	16
N. Protocol for purification of viral DNA (e.g. CMV) from urine	16
O. Protocol for purification of bacterial DNA (e.g. <i>Chlamydia trachomatis</i>) from cultures, biological fluids or clinical specimens	17
P. Protocol for purification of genomic DNA from insects.....	18
Q. Protocol for purification of genomic DNA from dental swabs	18
R. Protocol for purification of genomic DNA from buccal swabs.....	19
7. TROUBLESHOOTING.....	20
8. ORDERING INFORMATION.....	21
9. PRODUCT USE RESTRICTION AND WARRANTY	22

1. BASIC PRINCIPLE

SPEEDTOOLS TISSUE DNA EXTRACTION KIT is designed for the rapid isolation of highly pure genomic DNA from: **tissue, cells** (e.g. **bacteria, yeast**), **paraffin embedded samples, urine, feces, dried blood spots** and many other sources.

Lysis is achieved by incubation of the sample material in a proteinase K / SDS solution. Afterwards the lysate is treated with chaotropic ions and molecular grade ethanol to create the appropriate conditions for binding of DNA to the silica membrane of the column. The binding process is reversible and specific to nucleic acids. The washing steps efficiently remove contaminations and pure genomic DNA is finally eluted under low ionic strength conditions in a slightly alkaline elution buffer.

2. KIT CONTENTS

SPEEDTOOLS TISSUE DNA EXTRACTION KIT			
REAGENTS	10 Preps Ref. 21.135M	50 Preps Ref. 21.136	250 Preps Ref. 21.137
Buffer BT1	5 ml	20 ml	5 x 20 ml
Buffer BB3	3 ml	15 ml	5 x 15 ml
Buffer BB5 (concentrate)	3 ml	12 ml	5 x 12 ml
Buffer BBW	6 ml	30 ml	5 x 30 ml
Buffer BBE	3 ml	13 ml	5 x 13 ml
Proteinase K (lyophilized)	6 mg	30 mg	5 x 30 mg
Proteinase Buffer (PB)	400 µl	1.8 ml	5 x 1.8 ml
Speedtools Tissue Columns (plus Collection Tubes)	1 x 10	50	5 x 50
Collection Tubes	1 x 20	100	5 x 100

3. INTENDED USE

With the **SPEEDTOOLS TISSUE DNA** method, total DNA (genomic and mitochondrial) is prepared from tissue, cells and many other sources. It is also suitable for purify bacterial and viral DNA from these samples.

One Speedtools Tissue Column is capable of binding up to 60 µg of genomic DNA. It allows isolation of up to 35 µg of highly pure genomic DNA with an A260/280-ratio between 1.60 and 1.90.

For lysis of certain bacterial and yeast strains additional enzymes not included in the kit may be necessary. Relevant support protocols are provided with this manual.

The obtained DNA is ready to use for subsequent reactions like PCR, Southern blotting, or any kind of enzymatic reactions.

Table 1: General Characteristics of the Kit	
Sample material	□ 25 mg tissue/ 10 ² -10 ⁷ cultured cells
Expected yield	20-35 µg
Elution volume	60-100 µL
Binding capacity	60 µg
Time/ prep	20 min/prep (after lysis step)
Spin column	mini

4. PREPARATION OF WORKING SOLUTIONS AND STORAGE CONDITIONS

NOTE: Buffers BB3 and BBW contain guanidine hydrochloride therefore wear gloves and goggles.

All kit components can be stored at room temperature (18-25°C) through the expiration date printed on the packaging label.

Before starting any protocol with **SPEEDTOOLS TISSUE DNA EXTRACTION KIT** prepare the following reagents:

I. **Proteinase K:**

- ✓ **10 preps format:** Add **260 µL** of Proteinase Buffer to dissolve lyophilized Proteinase K.
- ✓ **50 preps format:** Add **1.35 mL** of Proteinase Buffer to dissolve lyophilized Proteinase K.

The Proteinase K Solution is stable at -20°C for at least 6 months.

II. **Buffer BB5:**

- ✓ **10 preps format:** Add **12 mL of molecular grade ethanol** (96-100%) to buffer BB5 concentrate. Mark the label of the bottle to indicate that ethanol was added.
- ✓ **50 preps format:** Add **48 mL of molecular grade ethanol** (96-100%) to buffer BB5 concentrate. Mark the label of the bottle to indicate that ethanol was added.

Wash Buffer BB5 can be stored at room temperature (18-25°C) for at least 12 months.

Upon storage, especially at low temperatures, a white precipitate may form in Buffer BB3. Dissolve such precipitates by incubation of the bottle at 50-70°C before use.

NOTE: During storage, especially at low temperatures, a white precipitate may form in Buffers BT1 or Buffer BB3. Such precipitates can be easily dissolved by incubating the bottle at 50-70 °C before use.

5. DNA ELUTION PROCEDURES

The *Standard Protocol* described in Section 6 has a recovery rate about 70-90%. It is possible to adapt the elution method and volume of elution buffer to the subsequent application of interest. In addition to the standard elution procedure there are several modifications possible. Use elution Buffer BBE for one of the following procedures:

- **High yield:** Perform two elution steps with the volume indicated in the individual protocol. About 90-100% of bound nucleic acid can be eluted.
- **High concentration:** Perform one elution step with 60% of the volume indicated in the individual protocol. Concentration of DNA will be approximately 30% higher than with standard elution. The yield of eluted DNA will be about 80%.
- **High yield and high concentration:** Apply half the volume of elution buffer indicated in the individual protocol, incubate for 3 min and centrifuge. Apply a second aliquot of elution buffer, incubate and centrifuge again. About 85-100% of bound nucleic acid is eluted in the standard elution volume at a high concentration.
- **Elution at 70 °C:** For certain samples heating the elution buffer to 70 °C increases the DNA yield.

Elution may also be performed with Tris-EDTA buffer (TE) of pH \geq 8. This will increase DNA stability especially during long term and/or multi use storage at 4°C or ambient temperature by inhibition of omnipresent DNases. However, EDTA interferes (depending on the final concentration) with certain downstream applications.

NOTE: Elution Buffer BBE (5 mM Tris/HCl, pH 8.5) provided with the kit does not contain EDTA.




For optimal performance of isolated DNA in downstream applications we recommend eluting with the supplied elution Buffer BBE and storage, especially long term, at -20°C. Several freeze-thaw cycles will not interfere with most downstream applications. Possible exceptions are detection of trace amounts of DNA or long-range PCR (e.g. > 10kb). Multiple freeze-thaw cycles or prolonged storage of eluted DNA at 4°C or room temperature may influence detection sensitivities or reaction efficiencies due to DNA shearing or adsorption to surfaces.

6. INSTRUCTION FOR USE





A. STANDARD PROTOCOL (for human or animal tissue and cultured cells)

Before starting the preparation:

- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Set an incubator or water bath to 56°C and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>Tissue Cut 25 mg human or animal tissue into small pieces. Place the sample in a microcentrifuge tube. Proceed with Step 2.</p> <p>Samples that are difficult to lyse can be ground under liquid nitrogen or may be treated in a mechanical homogenizer (Polytron, Ultra Turrax): Add 25 mg of tissue to a 1.5 ml centrifuge tube, add 50–75 µL phosphate buffered saline (PBS¹) and homogenize.</p> <p>Cultured cells Resuspend up to 10⁷ cells in a final volume of 200 µL Buffer BT1. Add 25 µL Proteinase K solution and 200 µL Buffer BB3. Vortex to mix and incubate the sample at 70°C for 10-15 min. Proceed with Step 4.</p>		PREPARE SAMPLE
2	<p>PRE-LYSIS</p> <p>Add 180 µL Buffer BT1 and 25 µL Proteinase K solution. Vortex to mix. Be sure that the samples are completely covered with lysis solution.</p> <p><i>If processing several samples, Proteinase K and Buffer BT1 may be premixed directly before use. Never mix Buffer BT1 and Proteinase K more than 10–15 min before addition to the sample: Proteinase K tends to self-digestion in Buffer BT1 without substrate.</i></p> <p>Incubate at 56°C until complete lysis is obtained (at least 1–3 h). Vortex occasionally during incubation or use a shaking incubator.</p> <p><i>Samples can be incubated overnight as well. If RNA-free DNA is crucial for downstream applications, a RNase digest may be performed: Add 20 µL RNase A (10 mg/ml) solution (not included) and incubate for an additional 5 min at room temperature.</i></p>		<p>+ 180 µL BUFFER BT1 + 25 µL PROTEINASE K</p> <p>Incubate 56°C 1-3h</p> <p>Or</p> <p>56°C overnight</p> <p>Vortex</p>
3	<p>SAMPLE LYSIS</p> <p>Vortex the samples. Add 200 µL of lysis Buffer BB3 to the samples and vortex the mixture vigorously. Incubate samples at 70°C for 10 min. Vortex briefly.</p> <p><i>If insoluble particles are visible, centrifuge for 5 min at high speed (e.g. 11,000 x g) and transfer the supernatant to a new microcentrifuge tube.</i></p>		<p>200 µL BUFFER BB3 + Vortex Incubate 70°C, 10 min</p>
4	<p>ADJUST DNA BINDING CONDITIONS</p> <p>Add 210 µL of molecular grade ethanol (96-100%) to each sample and vortex vigorously.</p> <p><i>After addition of ethanol a stringy precipitated may become visible. This will not affect the DNA isolation. Be sure to load all of the precipitate on the column in the following step.</i></p>		<p>+ 210 µL MOLECULAR GRADE ETHANOL Vortex</p>



¹ Buffer PBS sterile: dissolve 8 g NaCl; 0.2 g KCl, 1.44 g Na₂HPO₄; 0.24 g KH₂PO₄ in 800 ml H₂O. Adjust pH to 7.4 with HCl. Add H₂O to 1 Liter. Sterilise the Buffer PBS in the autoclave.

5	<p>BIND DNA</p> <p>For each sample, place one column into a 2 ml collection tube and apply the sample to the column. Centrifuge at 11,000 × g for 1 min. Discard the flow-through and place the column back into the collection tube.</p> <p><i>If the samples are not drawn completely through the matrix, repeat the centrifugation step at 11,000 × g. Discard flow-through.</i></p>		<p>Load lysate into a column</p> <p>1 min, 11,000 × g</p>
6	<p>WASH SILICA MEMBRANE</p> <ul style="list-style-type: none"> • Wash 1 <p>Add 500 µL Buffer BBW. Centrifuge 1 min at 11,000 x g. Discard the flow-through and place the column into the collection tube.</p> <ul style="list-style-type: none"> • Wash 2 <p>Add 600 µL Buffer BB5 to the column. Centrifuge 1 min at 11,000 x g. Discard the flow-through and place the column back into the collection tube.</p>		<p>+ 500 µL BUFFER BBW</p> <p>1 min, 11,000 × g</p> <p>+ 600 µL BUFFER BB5</p> <p>1 min, 11,000 × g</p>
7	<p>DRY SILICA MEMBRANE</p> <p>Place the column into 1.5 ml microcentrifuge tube and centrifuge 1 min at 11,000 x g.</p> <p><i>Residual ethanol is removed during this step.</i></p>		<p>1 min, 11,000 × g</p>
8	<p>ELUTE HIGHLY PURE DNA</p> <p>Place the column into a new 1.5 ml microcentrifuge tube and add 100 µL Buffer BBE. Dispense buffer directly onto the silica membrane. Incubate at room temperature for 1 min. Centrifuge 1 min at 11,000 x g.</p> <p><i>For alternative DNA elution procedures see section 5.</i></p>		<p>+ 100 µL BUFFER BBE (70°C)</p> <p>Incubate RT 1 min</p> <p>1 min, 11,000 × g</p>

B. Protocol for mouse or rat tails

Before starting the preparation:


- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Set an incubator or water bath to 56°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>Cut two 0.6 cm-pieces of mouse tail or one 0.6 cm-pieces of rat tail and place them in a 1.5 ml tube.</p>		CUT 2 PIECES (0.6cm)
2	<p>PRE-LYSIS</p> <p>Add 180 µL Buffer BT1 and 25 µL Proteinase K. Vortex to mix and incubate at 56°C overnight or until complete lysis is obtained.</p> <p><i>Lysis time can substantially be reduced down to approximately one hour if the tissue is broken up mechanically (e.g. if the tissue is cut into very small pieces before lysis).</i></p> <p>Vortex occasionally during incubation or use a shaking incubator. To remove residual bones or hair, centrifuge for 5 min at high speed (e.g. 11,000 x g). Transfer 200 µL supernatant to a new tube.</p> <p><i>If processing several samples, Proteinase K and Buffer BT1 may be premixed directly before use. Never mix Buffer BT1 and Proteinase K more than 10–15 min before addition to the sample: Proteinase K tends to self-digestion in Buffer BT1 without substrate.</i></p>		<p>+ 180 µL BUFFER BT1 + 25 µL PROTEINASE K</p> <p>Incubate 56°C, overnight Vortex</p> <p>Transfer 200 µL</p>
3	<p>SAMPLE LYSIS</p> <p>Add 200 µL of lysis Buffer BB3 to the samples and vortex the mixture vigorously.</p> <p><i>Buffer BB3 and molecular grade ethanol can be premixed before addition to the lysate.</i></p>		<p>+ 200 µL BUFFER BB3 Vortex</p>
4	<p>Proceed with Step 4 of the Standard Protocol.</p>		

C. Protocol for bacteria

Before starting the preparation:




- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Set an incubator or water bath to 56°C and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>Centrifuge ≤ 1 ml culture for 5 min at 8,000 x g. Remove supernatant.</p> <p><i>Depending on e.g. density of culture, culture medium, and bacterial strain, up to 1 ml of bacterial culture should be used.</i></p>		BACTERIA CULTURE (1ml maximum)
2	<p>PRE-LYSIS</p> <p>Resuspend the pellet in 180 μL Buffer BT1 by pipetting up and down. Add 25 μL Proteinase K. Vortex vigorously and incubate at 56°C until complete lysis is obtained (at least 1-3 h). Vortex occasionally during incubation or use a shaking incubator. Samples can be incubated overnight.</p> <p><i>If RNA-free DNA is crucial for downstream applications, an RNase digestion may be performed: Add 20 μL RNase A (10 mg/ml) solution (not included) and incubate for an additional 5 min at room temperature.</i></p> <p><i>Some strains, in particular Gram-positive bacteria, are more difficult to lyse. In such cases, a pre-incubation with a lytic enzyme is necessary. Resuspend the pelleted cells in 20 mM Tris/HCl; 2 mM EDTA; 1% Triton X-100; pH 8 (instead of Buffer BT1) supplemented with 20 mg/ml lysozyme or 0.2 mg/ml lysostaphin and incubate for 30-60 min at 37°C. Add 25 μL Proteinase K, incubate at 56 °C until complete lysis is obtained.</i></p>		<p>+ 180 μL BUFFER BT1 + 25 μL PROTEINASE K</p> <p>Vortex Incubate 56°C, 1-3h/overnight</p>
3	Proceed with Step 3 of the Standard Protocol.		

D. Protocol for yeast

Before starting the preparation:

- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Check that 10 mM EDTA, pH 8.0 is available.
- Check that sorbitol buffer² and lyticase or zymolase³ (not provided with the kit) is available for sample pre-lysis.
- Set an incubator or water bath to 30 °C; 56°C and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>Harvest 3 ml YPD yeast culture ($OD_{600} \leq 10$) by centrifugation for 10 min at 5,000 x g. Wash the cells once with 1 ml 10 mM EDTA, pH 8. Remove the supernatant and pellet the cells by centrifugation (5,000 x g, 10 min).</p>		<p>YEAST CULTURE (3ml YPD)</p> <p>10 min, 5,000×g</p> <p>Wash the cells in 1 ml 10mMEDTA,pH8</p>
2	<p>PRE LYSIS</p> <p>Resuspend the pellet in 600 µL sorbitol buffer. Add 50 U lyticase or zymolase. Incubate at 30°C for 30 min. This step degrades the yeast cell wall creating spheroplasts (formation of spheroplast may be checked microscopically).</p> <p>Centrifuge the mixture for 10 min at 2,000 x g, remove supernatant and resuspend the pelleted spheroplasts in 180 µL Buffer BT1. Add 25 µL Proteinase K and vortex vigorously.</p> <p>Incubate at 56°C until complete lysis is obtained (at least 1-3 h). Vortex occasionally during incubation or use a shaking incubator. Samples can be incubated overnight as well.</p> <p><i>If RNA-free DNA is crucial for downstream applications, an RNase digestion may be performed: Add 20 µL RNase A (20 mg/ml) solution (not included) and incubate for an additional 5 min at room temperature.</i></p>	 	<p>PELLET</p> <p>600 µL sorbitol buffer + 50 U lyticase/zymolase Incubate 30°C, 30 min</p> <p>10 min, 2,000×g + 180 µL BUFFER BT1 + 25 µL PROTEINASE K</p> <p>Incubate 56°C, 1-3 h/overnight Vortex</p>
3	Proceed with Step 3 of the Standard Protocol.		


² Buffer sorbitol: 1.2 M sorbitol; 10 mM CaCl₂; 0.1 M Tris/HCl pH 7.5; 35 mM β-Mercaptoetanol.

³ Other protocols use 5-200 U lyticase or zymolase depending on enzyme quality or brand. Increasing the enzyme concentration may be required if spheroplasts are not formed.

E. Protocol for dried blood spots (e.g. Guthrie cards; NucleoCards; and FTA® cards)

Before starting the preparation:




- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Set an incubator or water bath to 94°C and 56°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>Cut one or two dried blood spots as accurately as possible. Cut spots into small pieces and place them in a 1.5 ml centrifuge tube.</p> <p><i>The area of the dried blood spots should be between 15 and 30 mm².</i></p>		SAMPLE
2	<p>PRE-LYSIS</p> <p>Add 180 µL Buffer BT1 and mix by vortexing. Spin the samples briefly and place them in a water bath or heating block and heat 10 min at 94°C.</p> <p>Let the sample cool down. Add 25 µL Proteinase K, vortex and incubate at 56°C for 1h. Vortex occasionally during incubation or use a shaking water bath.</p> <p><i>Be sure that the samples are completely covered with lysis buffer during incubation.</i></p>		<p>+ 180 µL BUFFER BT1 Incubate 94°C, 10 min</p> <p>Cool + 25 µL PROTEINASE K Incubate 56°C, 1 h Vortex</p>
3	<p>SAMPLE LYSIS</p> <p>Add 200 µL Buffer BB3 vortex vigorously and incubate at 56°C for 10 min.</p>		<p>200 µL BUFFER BB3 Vortex Incubate 56°C, 10 min</p>
4	<i>Proceed with Step 4 of the Standard Protocol.</i>		

F. Protocol for genomic DNA and viral DNA from blood samples

Before starting the preparation:


- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Set an incubator or water bath to 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1-2	PREPARE SAMPLE AND PRE-LYSIS These steps are not necessary for the present protocol.		Not necessary
3	SAMPLE LYSIS Add 25 µL Proteinase K and up to 200 µL blood , buffy coat, or body fluid sample (equilibrated to room temperature). <i>For sample volumes less than 200 µL, add PBS to adjust the volume to 200 µL. If purifying DNA viruses, we recommend starting with 200 µL serum or plasma. If cultured cells are used, resuspend up to 5 x 10⁶ cells in a final volume of 200 µL PBS.</i> Add 200 µL Buffer BB3 to the samples and vortex the mixture vigorously (10-20 sec). Incubate samples at room temperature for 5 min . Mix and incubate sample at 70 °C for 10-15 min . <i>The lysate should become brownish during incubation with Buffer BB3. Increase incubation time with Proteinase K (up to 30 min) and vortex one or twice vigorously during incubation if processing older or clotted blood samples.</i>		25 µL PROTEINASE K + 200 µL sample + 200 µL BUFFER BB3 Incubate RT, 5 min Incubate 70°C, 10-15 min
4	ADJUST DNA BINDING CONDITIONS Add 210 µL molecular grade ethanol (96-100 %), to each sample, and vortex again.		210 µL Molecular Grade Ethanol Vortex
5	BIND DNA For each sample, place one column into a 2 ml collection tube and apply the sample to the column. Centrifuge at 11,000 × g for 1 min . Discard the flow-through and place the column back into the collection tube. <i>If the samples are not drawn completely through the matrix, repeat the centrifugation step at 11,000 × g. Discard flow-through.</i>		Load lysate into a column 1 min, 11,000 × g
6	Proceed with Step 4 of the Standard Protocol.		

G. Protocol for hair roots

Before starting the preparation:




- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Set an incubator or water bath to 56°C and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	PREPARE SAMPLE Cut off the hair roots from the hair sample (up to 100) and collect them in 1.5 ml collection tube.		HAIR ROOTS
2	PRE-LYSIS Add 180 µL Buffer BT1 to the hair roots and freeze the samples in liquid nitrogen. Thaw samples in a 56°C water bath. Repeat this procedure 4 times. Add 25 µL Proteinase K , mix by vortexing, and incubate 6-8 h or overnight at 56°C . Use a shaking water bath/incubator or vortex occasionally.		+ 180 µL BUFFER BT1 5 cycles of freeze/thaw with liquid Nitrogen 25 µL PROTEINASE K Incubate 56°C, 6-8 h/overnight Vortex
3	<i>Proceed with Step 3 of the Standard Protocol.</i>		

H. Protocol from paraffin-embedded tissue

Before starting the preparation:


- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Check that n-octane or xylene is available.
- Set an incubator or water bath to 37°C; 56°C; and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	PREPARE SAMPLE Prepare small sections (up to 25 mg) from blocks of fixed, embedded tissue . If possible, trim excess paraffin from the block before slicing. Handle the sections with tweezers or toothpicks and place the samples into microcentrifuge tubes. Add 1 ml n-octane or xylene to each tube. Vortex vigorously and incubate at room temperature for about 30 min . Vortex occasionally. Centrifuge at 11,000 x g for 3 min . Pipette off supernatant. Add 1 mL molecular grade ethanol (96-100%) to each tube. Close and mix by inverting several times. Centrifuge at 11,000 x g for 3 min . Pipette off supernatant. Repeat the ethanol washing step. Pipette off as much of the ethanol as possible. Incubate the open tube at 37°C until the ethanol has evaporated (~15 min).	  	SMALL SECTIONS + 1 mL n-OCTANE Vortex Incubate RT, 30 min 3 min, 11,000×g + 1 mL ETANOL 3 min, 11,000×g + 1 mL MOLECULAR GRADE ETHANOL 3 min, 11,000×g Incubate open tube at 37°C
2	<i>Proceed with Step 2 of the Standard Protocol.</i>		

I. Protocol for detection of *Mycobacterium tuberculosis* or *Legionella pneumophila* in sputum or bronchoalveolar lavage

Before starting the preparation:

- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Prepare N-acetyl cystein/NaOH⁴.
- Set an incubator or water bath to 56°C and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.


STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>Add 200-500 µL sputum or bronchoalveolar lavage to an equal volume N-acetyl cystein/NaOH. Vortex gently to mix.</p> <p>Incubate the mixture for 25 min at room temperature with shaking.</p> <p>Adjust the volume to 25 mL with sterile water. Centrifuge for 30 min at 4,000 x g. Discard the supernatant.</p> <p>Resuspend the pellet in 0.5-1 mL Buffer BT1, depending on sample viscosity (Kit only includes 20 mL of Buffer BT1).</p> <p>Transfer 200 µL of the resuspended sample to a new microcentrifuge tube.</p> <p><i>For this protocol depending on the number and the viscosity of the samples you may require and additional bottle of Buffer BT1 (BIOTOOLS Ref. 21.161).</i></p>		<p>ESPUTUM 200-500 µL + equal vol N-acetyl cystein/NaOH</p> <p>Vortex Incubate RT 25 min</p> <p>Adjust vol to 25 mL with sterile H₂O 30 min, 4,000 x g</p> <p>Resuspend pellet in 0.5-1 mL BUFFER BT1</p> <p>Transfer 200 µL</p>
2	Proceed with Step 2 of the Standard Protocol.		

J. Protocol for extraction of genomic DNA from stool

For this protocol depending on the number and the quality of the samples you may require and additional bottle of Buffer BT1 (BIOTOOLS Ref. 21.161).

Before starting the preparation:

- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Check if TE buffer⁵ is available.
- Set an incubator or water bath to 56°C and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>Add 250 mg feces to 1 mL TE buffer. Resuspend the sample by vigorously vortexing (30 s).</p> <p>Centrifuge the sample for 15 min at 4,000 x g. Discard the supernatant. Resuspend the pellet in 0.2-1 mL Buffer BT1. Use as much buffer as necessary for good resuspension of the sample.</p> <p><i>The prepared pellet contains, among other constituents, cells from the digestive tract and bacteria.</i></p> <p>Transfer 200 µL of the resuspended sample to a new microcentrifuge tube.</p>		<p>FECES (250 mg) + 1 mL buffer TE</p> <p>Vortex</p> <p>15 min, 4,000 x g</p> <p>Resuspend in 0.2-1 mL BUFFER BT1</p> <p>Transfer 200 µL</p>
2	Proceed with the addition of 25 µL Proteinase K in Step 2 of the Standard Protocol.		

⁴ N-acetyl cystein/NaOH: 2 g NaOH; 1.45 g sodium citrate; 0.5 g N-acetyl cystein. Add water to 100 ml.







⁵ Buffer TE: 10 mM Tris/HCl; 1 mM EDTA, pH 8.



K. Protocol for viral DNA from stool

For this protocol you will require additional bottles of Proteinase K solution (BIOTOOLS Ref. 20050).

Before starting the preparation:

- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Prepare 0.9 % NaCl.
- Set an incubator or water bath to 70 °C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>Suspend the stool sample in 0.9 % NaCl.</p> <p>Centrifuge aliquots of the stool sample (1 mL) for 5 min at 800 × g at room temperature.</p> <p>Carefully reunite supernatant (do not touch the pellet).</p> <p>Filtrate supernatant through 0.22-0.45 µm.</p> <p>Fractionate the filtrate and centrifuge for 1 min at 11,000 × g.</p>		<p>PREPARE SAMPLE + 0.9% NaCl</p> <p>1 min, 800 × g</p> <p>Filtrate sample</p> <p>1 min, 11,000 × g</p>
2	<p>PRE-LYSIS</p> <p>Carefully remove the supernatant by decanting. Add 400 µL Buffer BT1 and 35 µL Proteinase K solution. Vortex to mix.</p>		<p>+ 400 µL BUFFER BT1 + 35 µL PROTEINASE K</p>
3	<p>SAMPLE LYSIS</p> <p>Vortex the samples. Add 400 µL of lysis Buffer BB3 to the samples and vortex the mixture vigorously. Incubate samples at 70°C for at least 30 min. Vortex briefly.</p>		<p>400 µL BUFFER BB3 + Vortex Incubate 70°C, 30 min</p>
4	<p>ADJUST DNA BINDING CONDITIONS</p> <p>Add 420 µL of molecular grade ethanol (96-100%) to each sample and vortex vigorously.</p>		<p>+ 420 µL MOLECULAR GRADE ETHANOL Vortex</p>
5	<p>BIND DNA</p> <p>For each sample, place one column into a 2 ml collection tube and apply the sample to the column. Centrifuge at 4,500 × g for 1 min. Discard the flow-through and place the column back into the collection tube.</p> <p><i>If the samples are not drawn completely through the matrix, repeat the centrifugation step at 11,000 × g. Discard flow-through.</i></p>		<p>Load lysate into a column</p> <p>1 min, 4,500 × g</p>
6	<p>WASH SILICA MEMBRANE</p> <ul style="list-style-type: none"> • Wash 1 <p>Add 600 µL Buffer BBW. Centrifuge 1 min at 4,500 x g. Discard the flow-through and place the column into the collection tube.</p> <ul style="list-style-type: none"> • Wash 2 <p>Add 600 µL Buffer BB5 to the column. Centrifuge 1 min at 4,500 x g. Discard the flow-through and place the column back into the collection tube.</p> <ul style="list-style-type: none"> • Wash 3 <p>Add 600 µL Buffer BB5 to the column. Centrifuge 2 min at 11,000 x g. Discard the flow-through and place the column back into the collection tube.</p>		<p>+ 600 µL BUFFER BBW 1 min, 4,500 × g</p> <p>+ 600 µL BUFFER BB5 1 min, 4,500 × g</p> <p>+ 600 µL BUFFER BB5 2 min, 11,000 × g</p>




7	<p>DRY SILICA MEMBRANE</p> <p>Place the column into 1.5 ml microcentrifuge tube and incubate with open lid for 1-2 min at 70 °C.</p> <p><i>Residual ethanol is removed during this step.</i></p>		<p>Incubate 70°C, 1-2 min</p>
8	<p>ELUTE HIGHLY PURE DNA</p> <p>Place the column into a new 1.5 ml microcentrifuge tube and add 100 µL Buffer BBE. Dispense buffer directly onto the silica membrane. Incubate with closed lid for 3-5 min.</p> <p>Centrifuge 1 min at 4,500 x g.</p> <p><i>For alternative DNA elution procedures see section 5.</i></p>		<p>+ 100 µL BUFFER BBE</p> <p>Incubate RT 3-5 min</p> <p>1 min 4,500 x g</p>

L. Protocol for detection of EHEC bacteria in food (e.g. fresh cows' milk)

In humans, Vero toxin-forming *E. coli* strains (VTEC, EHEC) can cause diseases. The main reservoirs and sources of infection for humans are horned cattle and the corresponding foods, especially raw or insufficiently cooked minced meat and raw milk.

Before starting the preparation:

- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Prepare mTSB⁶ and 3.2 M sodium acetate.
- Check if Novobiocin is available.
- Set an incubator or water bath to 37 °C; 56°C and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.


STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>To a sterile 1 liter flask, add 25 mL milk and 225 mL pre-warmed (37°C) mTSB medium (supplied with Novobiocin). Incubate the mixture in a shaking water bath for 5-6 h or overnight at 37°C.</p> <p>Centrifuge 100 mL culture for 40 min at 6,000 x g.</p> <p>Gently pour off the supernatant and resuspend the pellet in 2 mL sterile water.</p> <p>Centrifuge for 10 min at 10,000 x g.</p>		<p>MILK (25 mL) + 225 mL sterile mTSB</p> <p>Incubate 37°C with shaking, 5-6 h/overnight</p> <p>100 mL culture 40 min, 6,000 x g</p> <p>Resuspend cells in 2 mL sterile H₂O</p> <p>10 min, 10,000 x g</p>
2	<p>PRE-LYSIS</p> <p>Resuspend the pellet in 180 µL Buffer BT1 and add 25 µL Proteinase K.</p>		<p>+ 180 µL BUFFER BT1</p> <p>+ 25 µL PROTEINASE K</p>
3-8	<p><i>Proceed with Step 3 of the Standard Protocol.</i></p> <p><i>Perform two elution steps with the indicated elution volume of the Standard Protocol.</i></p>		200 µL ELUATE
9	<p>After elution of the DNA, proceed with the following step.</p> <p>Precipitate the obtained DNA by adding 20 µL of 3.2 M sodium acetate and 400 µL of molecular grade ethanol to 200 µL the eluate.</p> <p>Centrifuge for 30 min at 11,000 x g. Discard supernatant and wash the pellet with 1 mL 70% molecular grade ethanol and resuspend in 10 µL sterile water.</p>		<p>+ 20 µL 3.2 M SODIUM ACETATE</p> <p>+ 400 µL MOLECULAR GRADE ETHANOL</p> <p>30 min, 11,000 x g</p> <p>Wash pellet with 1 mL MOLECULAR GRADE ETHANOL 70%</p> <p>Resuspend in 10 µL sterile H₂O</p>

⁶ mTSB Medium: 30 g Tryptic soy broth (Gibco); 1.5 g biles salts No. 3 (Oxoid); 1.5 g KH₂PO₄. Add 900 ml de H₂O. Filter the medium and adjust the pH with 2 M NaOH to 7.4. Add water to 1 liter and autoclave for 15 min at 121°C.

M. Protocol for purification of bacterial DNA (e.g. *Borrelia burgdorferi*) from urine

Before starting the preparation:




- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Set an incubator or water bath to 56°C and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.





STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>Centrifuge 1 mL urine sample at 13,000 x g for 30 min. Discard the supernatant, add again 1 mL urine sample to the pellet and centrifuge at 13,000 x g for 30 min. Repeat this step a third time.</p> <p><i>The sample material should be fresh and storage at -20 to -80°C is only recommendable for a couple of days. After thawing incubate the sample at 40°C as long as all precipitates, are dissolved. Urine tends to form precipitates when stored at low temperatures. If a complete solution does not happen let the precipitate sediment and proceed with step 1 of the protocol using only the supernatant</i></p>		<p>1 mL URINE 30 min, 13,000 x g</p> <p>+</p> <p>1 mL URINE 30 min, 13,000 x g</p> <p>+</p> <p>1 mL URINE 30 min, 13,000 x g</p>
2	Proceed with the Step 2 of the Standard Protocol.		

N. Protocol for purification of viral DNA (e.g. CMV) from urine

Before starting the preparation:

- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Set an incubator or water bath to 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.


STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>Prepare 3-4 aliquots (1 mL each aliquot) of urine sample. Centrifuge aliquots for 10 min at 13,000 x g. Carefully decant the supernatant of the aliquots.</p> <p>If frozen urine samples are used precipitates may appear after defrosting, which must be dissolved before the centrifugation step. This can be done through a 30 min incubation step at 37-40°C. If a complete solution does not happen let the precipitate sediment and proceed with Step 1 of the protocol using only the supernatant.</p>		<p>SAMPLE divided in aliquots of 1 mL</p> <p>1 mL URINE 10 min, 13,000 x g</p>
2	<p>PRE-LYSIS</p> <p>Resuspend the pellet in 180 µL Buffer BT1 and 25 µL Proteinase K.</p> <p>Resuspend the first pellet with reagents indicated above. Transfer the resuspended solution of the first tube to the second tube and the resuspended solution of the second tube to the third tube and so on. Finally continue with step 3.</p>		<p>+</p> <p>180 µL BUFFER BT1</p> <p>+</p> <p>25 µL PROTEINASE K</p>
3	<p>SAMPLE LYSIS</p> <p>Add 200 µL of lysis Buffer BB3 and incubate at least 20 min at 70°C. Vortex briefly.</p>		<p>200 µL BUFFER BB3</p> <p>Vortex and Incubate 70°C, 20 min</p>

4	<p>ADJUST DNA BINDING CONDITIONS</p> <p>Add 210 µL of molecular grade ethanol (96-100%) to the sample and vortex vigorously.</p>		<p>+ 210 µL MOLECULAR GRADE ETHANOL Vortex</p>
5	<p>BIND DNA</p> <p>For each sample, place one column into a 2 mL collection tube and apply the sample to the column. Centrifuge at 4,500 × g during 1 min. Discard the flow-through and place the column back into the collection tube.</p>		<p>Load lysate into a column 1 min, 4,500 × g</p>
6	<p>WASH SILICA MEMBRANE</p> <ul style="list-style-type: none"> • 1st Wash <p>Add 500 µL Buffer BBW. Centrifuge 1 min at 4,500 × g. Discard flow-through and place the column back into the collection tube.</p> <ul style="list-style-type: none"> • 2nd Wash <p>Add 600 µL Buffer BB5. Centrifuge 2 min at 11,000 × g. Discard the flow-through and place the column back into the collection tube.</p>		<p>+ 500 µL BUFFER BBW 1 min, 4,500 × g</p> <p>+ 600 µL BUFFER BBW 2 min, 11,000 × g</p>
7	<p>DRY SILICA MEMBRANE</p> <p>Incubate with open lid for 1-2 min at 70°C. <i>Residual ethanol is removed during this step.</i></p>		<p><i>Incubate tubes with open lid Incubate 70° 1-2 min</i></p>
8	<p>ELUTE HIGHLY PURE DNA</p> <p>Add 70 µL Buffer BBE, close the lid and incubate for further 3-5 min. Centrifuge 1 min at 4,500 × g.</p>		<p>+ 70 µL BUFFER BBE Incubate 70°C 3-5 min 1 min, 4,500 × g</p>

O. Protocol for purification of bacterial DNA (e.g. *Chamydia trachomatis*) from cultures, biological fluids or clinical specimens

Before starting the preparation:

- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Set an incubator or water bath to 56°C and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <ul style="list-style-type: none"> • Isolation of bacterial DNA from bacterial cultures or biological fluids: pellet bacteria by centrifugation for 5 min at 13,000 × g. • Isolation of bacterial DNA from eye, nasal or pharyngeal swabs: collect samples, add 2 mL PBS⁷ containing a common fungicide, and incubate for several hours at room temperature. Pellet bacteria by centrifugation for 5 min at 13,000 × g 		<p>SAMPLE 5 min, 13,000 × g</p>
2	<i>Proceed with Step 2 of the Standard Protocol.</i>		

⁷ Buffer PBS sterile: dissolve 8 g NaCl; 0.2 g KCl, 1.44 g Na₂HPO₄; 0.24 g KH₂PO₄ in 800 ml H₂O. Adjust pH to 7.4 with HCl. Add H₂O to 1 Liter. Sterilise the Buffer PBS in the autoclave.

P. Protocol for purification of genomic DNA from insects

Before starting the preparation:


- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Set an incubator or water bath to 56°C and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	PREPARE SAMPLE Homogenize not more than 50 mg insects under liquid nitrogen and transfer the powder into a 1.5 mL centrifuge tube.		SAMPLE (<50mg) Homogenize under liquid nitrogen
2	<i>Proceed with the Step 2 of the Standard Protocol.</i>		

Q. Protocol for purification of genomic DNA from dental swabs

Before starting the preparation:




- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Set an incubator or water bath to 70 °C and 95 °C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	PREPARE SAMPLE Place swab material (paper, cotton, brushes, plastic) in a 1.5 mL microcentrifuge tube.		SAMPLE
2	PRE-LYSIS Add 180 µL Buffer BT1 and 25 µL Proteinase K to each sample. Close the microcentrifuge tube and spin briefly for 15 s at 1,500 x g in order to submerge the swab material completely. Incubate at room temperature for 5 min . Vortex the tube vigorously for 15 s and spin briefly for 15 s at 1,500 x g . Incubate the tubes at 70°C in an incubator for 10 min . Place a weight on top of the tube in order to prevent the caps from popping off. Shift the temperature to 95°C for 5 min . Spin briefly for 15 s at 1,500 x g to collect any sample from the lids. Open the microcentrifuge tubes. <i>Depending on the bacterial strains that are to be detected, incubation at 95°C can be skipped.</i>		180 µL BUFFER BT1 + 25 µL PROTEINASE K Spin Incubate 5 min at TA vortex and spin Incubate 10 min, 70°C Incubate 5 min, 95°C Spin
3	SAMPLE LYSIS Transfer as much as possible of the lysate solution to a 1.5 mL microcentrifuge tube. Discard swab and continue with recovered solution.		
	<i>Proceed with Step 3 of the Standard Protocol.</i>		

R. Protocol for purification of genomic DNA from buccal swabs

Before starting the preparation:

- Check if Buffer BB5 and Proteinase K were prepared according to Section 4.
- Check if PBS⁸ is available.
- Set an incubator or water bath to 56°C and 70°C.
- Before elution and if necessary preheat Elution Buffer BBE to 70°C.

STEP	DESCRIPTION		
1	<p>PREPARE SAMPLE</p> <p>Collect the samples with cotton; dacron® (Daigger); or C.E.P. swabs (Gibco BRL). Scrape firmly against the inside of each cheek several times and let the swabs air dry.</p> <p><i>The individual should not have consumed food or drink within 30 min before collection of the sample.</i></p>		SAMPLE
2	<p>PRE-LYSIS</p> <p>Place the dry swab material in 2 mL microcentrifuge tube. Add 400-600 µL PBS and 25 µL Proteinase K to the swabs.</p> <p><i>The volume of PBS is depending on the type of swab used: for cotton and dacron swabs 400 µL are sufficient; for C.E.P. swabs 600 µL are necessary.</i></p> <p>Mix by vortexing 2 x 5 sec and incubate 10 min at 56°C.</p>		<p>+ 400-600 µL buffer PBS + 25 µL PROTEINASE K</p> <p>Incubate 10 min, 56°C</p>
3	<p>SAMPLE LYSIS</p> <p>Transfer as much as possible of the lysate solution to a 1.5mL microcentrifuge tube.</p> <p>Discard swab and continue with recovered solution.</p> <p>Add one volume Buffer BB3 (400-600 µL, depending on the swab type/volume of PBS buffer used) and vortex vigorously. Incubate the samples at 70°C for 10 min.</p>		<p>1 vcl BUFFER BB3</p> <p>Incubate 10 min, 70°C</p>
4	<p>ADJUST DNA BINDING CONDITIONS</p> <p>Add one volume (200-300 µL, depending on the swab type) molecular grade ethanol (96-100%) to each sample and mix by vortexing.</p>		<p>+ 1 vcl MOLECULAR GRADE ETHANOL Vortex</p>
5	<p>BIND DNA</p> <p>Transfer 600 µL of the samples from the 2 mL microcentrifuge tubes into the columns. Centrifuge at 11,000 x g for 1 min.</p> <p><i>If the samples are not drawn through completely, repeat the centrifugation.</i></p> <p>Discard the flow-through and place the column back into the collection tube. Repeat Step 5 once or twice, depending on the lysis volume.</p>		<p>Load lysate into a column</p> <p>1 min, 11,000 × g</p>
6	<p>When all of the lysate has been applied to the columns, discard Collection Tube and place the column in a new Collection Tube.</p> <p>Proceed with Step 6 of the Standard Protocol.</p>		

⁸ Buffer PBS sterile: dissolve 8 g NaCl; 0.2 g KCl, 1.44 g Na₂HPO₄; 0.24 g KH₂PO₄ in 800 ml H₂O. Adjust pH to 7.4 with HCl. Add H₂O to 1 Liter. Sterilise the Buffer PBS in the autoclave.

7. TROUBLESHOOTING

Problem	Possible cause and suggestions
No yield or poor DNA yield	<p>Incomplete cell lysis</p> <ul style="list-style-type: none"> • Sample not thoroughly homogenized and mixed with Buffer BT1 / Proteinase K. The mixture has to be vortexed vigorously immediately after addition of Buffer BT1. • Decreased Proteinase K activity: store dissolved Proteinase K at –20°C for at least 6 months. <p>Reagents not applied properly</p> <ul style="list-style-type: none"> • Prepare Buffers BB3, BB5, and Proteinase K according to instructions (section 4). Add molecular grade ethanol to lysates before loading them on the columns. <p>Suboptimal elution of DNA from the column</p> <ul style="list-style-type: none"> • Preheat Buffer BBE to 70°C before elution. Apply Buffer BBE directly onto the center of the silica membrane. • Elution efficiencies decrease dramatically if elution is done with buffers with pH < 7.0. Use slightly alkaline elution buffers like Buffer BBE (pH 8.5). • Especially when expecting high yields from large amounts of material, we recommend elution with 200 µL Buffer BBE and incubation of the closed columns in an incubator at 70°C for 5 min before centrifugation.
Poor DNA quality	<p>Incomplete lysis</p> <ul style="list-style-type: none"> • Sample not thoroughly homogenized and mixed with Buffer BT1 / Proteinase K. The mixture has to be vortexed vigorously immediately after addition of Buffer BT1. • Decreased Proteinase K activity: store dissolved Proteinase K at –20°C for at least 6 months. <p>Reagents not applied properly</p> <ul style="list-style-type: none"> • Prepare Buffers BB3, BB5, and Proteinase K according to instructions (section 4). Add molecular grade ethanol to lysates before loading them on the columns. <p>RNA in sample</p> <ul style="list-style-type: none"> • If DNA free of RNA is desired, add RNase A solution, not supplied with the kit, before addition of lysis Buffer BB3 and incubate at 37°C for 5 min.
Columns clogged	<p>Too much sample material used</p> <ul style="list-style-type: none"> • Do not use more sample material than recommended (25 mg for most tissue types). If insoluble material like bones or hair remains in the lysate, spin down the debris and transfer the clear supernatant to a new microcentrifuge tube before proceeding with addition of Buffer BB3 and molecular grade ethanol. <p>Incomplete lysis</p> <ul style="list-style-type: none"> • Sample not thoroughly homogenized and mixed with Buffer BT1 / Proteinase K. The mixture has to be vortexed vigorously immediately after addition of Buffer BT1. • Decreased Proteinase K activity: store dissolved Proteinase K at –20°C for at least 6 months. <p>Reagents not applied properly</p> <ul style="list-style-type: none"> • Prepare Buffer BB5 and Proteinase K according to instructions (section 4). Add molecular grade ethanol to lysates before loading them on the columns.

Suboptimal performance of genomic DNA in enzymatic reactions	<p>Carryover of ethanol or salt</p> <ul style="list-style-type: none"> • Be certain to centrifuge ≥ 1 min at 11,000 x g in order to remove all of ethanolic Buffer BB5 before eluting the DNA. If, for any reason, the level of Buffer BB5 has reached the column outlet after drying, repeat the centrifugation. • Do not chill Buffer BB5 before use. Cold buffer will not remove salt effectively. Equilibrate Buffer BB5 to room temperature before use. <p>Contamination of DNA with inhibitory substances</p> <ul style="list-style-type: none"> • DNA eluted with Tris/EDTA buffer (TE). EDTA may inhibit enzymatic reactions. Repurify DNA and eluate in Buffer BBE. • If the A260/280 ratio of the eluate is ≤ 1.6, repeat the purification procedure: add 1 vol Buffer BB3 plus 1 vol molecular grade ethanol (96-100%) to the eluate. Load the mixture onto a column and proceed with step 5 of the Standard Protocol.
--	---

8. ORDERING INFORMATION

SPEEDTOOLS KIT	10 PREPS	50 PREPS	250 PREPS
SPEEDTOOLS DNA EXTRACTION KIT	Ref. 21.130M	Ref. 21.131	Ref. 21.132
SPEEDTOOLS TISSUE DNA EXTRACTION KIT	Ref. 21.135M	Ref. 21.136	Ref. 21.137
SPEEDTOOLS RNA VIRUS EXTRACTION KIT	Ref. 21.140M	Ref. 21.141	Ref. 21.142
SPEEDTOOLS FOOD DNA EXTRACTION KIT	Ref. 21.175M	Ref. 21.176	Ref. 21.177
SPEEDTOOLS PLANT DNA EXTRACTION KIT	Ref. 21.170M	Ref. 21.171	Ref. 21.172
SPEEDTOOLS TOTAL RNA EXTRACTION KIT	Ref. 21.210M	Ref. 21.211	Ref. 21.212
SPEEDTOOLS PCR CLEAN-UP KIT	Ref. 21.200M	Ref. 21.201	Ref. 21.202
SPEEDTOOLS PLASMID DNA PURIFICATION KIT	Ref. 21.220M	Ref. 21.221	Ref. 21.222

9. PRODUCT USE RESTRICTION AND WARRANTY

1. Product for research purposes and *in vitro* uses only.
2. No claim or representation is intended for its use to identify any specific microorganism or for clinical use (diagnostic, prognosis, therapeutic, or blood bank).
3. It is rather the responsibility of the user to verify the use of the kit for a specific application range as the performance characteristic of this kit has not been verified for a specific microorganism.
4. The kit includes documentation stating specifications and other technical information. Follow the instructions in order to obtain highly pure genomic DNA.
5. BIOTOOLS warrants to meet the stated specifications of the kit. Any product not fulfilling the specifications included in the product sheet will be replaced.
6. BIOTOOLS obligation and purchaser's rights under this warranty are limited to the replacement by BIOTOOLS of any product that is shown defective in fabrication, and that must be returned to BIOTOOLS, freight prepaid, or at BIOTOOLS' option, replacement of the purchasing price.
7. BIOTOOLS does not warrant against damages or defects arising in shipping and handling (transport insurance for customers excluded). Any complaint on damaged goods during transport must be directed to the handling or transport agent.
8. BIOTOOLS has no responsibility for damages, whether direct or indirect, incidental or consequential of improper or abnormal use of this product. Nor has any responsibility for defects in products or components not manufactured by BIOTOOLS, or against damages resulting from such non-BIOTOOLS components or products.
9. BIOTOOLS makes no other warranty of any kind whatsoever, and specifically disclaims and excludes all other warranties of any kind or nature whatsoever, directly or indirectly, express or implied, including, without limitation, as to the suitability, reproductively, durability, fitness for a particular purpose or use, merchantability, condition, or any other matter with respect to BIOTOOLS products.
10. The warranty provided herein and the data, specifications and descriptions of this kit appearing in BIOTOOLS published catalogues and product literature are BIOTOOLS sole representations concerning the product and warranty. No other statements or representations, written or oral, by BIOTOOLS employees, agent or representatives, except written statements signed by a duly authorized officer of BIOTOOLS are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.
11. Product claims are subject to change. Therefore please contact our Technical Service Team for the most up-to-date information. You may also contact your local distributor for general scientific information.
12. Applications mentioned in BIOTOOLS literature are provided for informational purposes only. BIOTOOLS does not warrant that all applications have been tested in BIOTOOLS laboratories using BIOTOOLS products. BIOTOOLS does not warrant the correctness of any of those applications. For more information contact our Technical Dept (technicalsupport@biotools.eu).

Manufactured by:

BIOTOOLS, Biotechnological & Medical Laboratories, S.A. has been evaluated and certified to accomplish ISO 9001 requirements for the following activities: Research and development of biotechnology products and manufacture of biotechnology and in vitro products.

Valle de Tobalina – 52 – Nave 39, 28021 Madrid – Spain



© 2008 BIOTOOLS, Biotechnological & Medical Laboratories, S.A. All rights reserved.