

Random Primer

Product Description:

The primer is a mixture of single-stranded random hexanucleotides with 5'- and 3'-hydroxyl ends. Random hexamer primers can be used for first-strand cDNA synthesis and cloning. The primer is supplied in lyophilized form 30 µg.

Form: Supplied in lyophilized form

Storage: 4°C or -20°C

Pack size: 30 µg ~ 17 nmol

Notes for Use:

Stored at 4°C or -20°C when lyophilised; once the product is resuspended in nuclease-free water store at **-20°C**. For maximum recovery prior removing the cap of the tube spin it briefly to bring down the content. Reconstitute with nuclease-free water (included vial) and vortex briefly. Avoid repeat freeze-thaw cycles. Aliquots are stable for at least 24 months at -20°C. Further dilutions can be made in assay buffer.

Order Information:

Ref.	Product	
20.401	Random Primers	30 µg
20.402	Oligo dT(15)	30 µg

Notice to buyers/users:

Some of the applications which may be performed with this product are covered by applicable patents in certain countries. The purchase of this product does not include or provide a license to perform patented applications. Users may be required to obtain a license depending on the country and/or application.

Random Primer

Descripción del Producto

Mezcla de hexanucleótidos de cadena sencilla con grupos hidroxilo en sus extremos 5' y 3'. Puede utilizarse para síntesis de ADNc y para clonar. El producto se encuentra liofilizado (30 µg).

Formato: liofilizado

Almacenamiento: 4°C o -20°C

Tamaño del Producto: 30 µg ~ 17 nmol

Notas de Uso

Almacenar a 4°C o a -20°C el vial liofilizado; una vez que el vial se ha resuspendido en agua libre de nucleasas almacenar a **-20°C**. Antes de abrir el vial centrifugue brevemente para así recuperar todo su contenido. Reconstituya el vial con agua libre de nucleasas (vial incluido en el pedido) y mezcle con el vortex brevemente. Evite ciclos de congelación y descongelación. Las alícuotas son estables por 2 años a -20°C. Las diluciones posteriores pueden realizarse en el buffer de trabajo.

Información de Pedidos

Ref.	Producto	
20.401	Random Primers	30 µg
20.402	Oligo dT(15)	30 µg

Nota a compradores/usuarios:

Algunas de las aplicaciones que pueden ser realizadas con este producto están protegidas por patentes aplicables en ciertos países. La compra de este producto no incluye una licencia para llevar a cabo aplicaciones patentadas. Los usuarios podrían tener que adquirir una licencia dependiendo del país y/o aplicación.